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Sascha Seifert &  
Dominik Suslik  
info@vitalityflossing.de

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### General notes:

- Consult a doctor before use and only carry out exercises under the care of a medical professional.
- Improper use can lead to serious injury!
- Avoid sharp objects as holes/tears can lead to the band ripping.
- Store the band rolled up in a case, away from direct sunlight or heat (heating).
- This product contains pure latex (100% natural latex) that can cause allergic reactions.
- Not a toy! Not suitable for children!
- Avoid contact with wounds and mucus membranes as well as long-term contact with skin.
- Clean with alcohol-based disinfectant cloths / water & soap

### Practical Use

#### Thumb:

#### Wrapping:

Fixation in the centre of the back of the hand.

Wrap around the surface of the hand and, finally, the thumb joint.

Secure the ends of the band at the wrist.

#### Active activation:

Move, bend and stretch the thumb actively in all directions.

#### Passive activation:

Supported thumb movements

#### Elbow:

#### Wrapping:

Fixation underneath the elbow.

Wrap over the elbow towards the shoulder.

Secure the ends of the band at the upper arm.

#### Activation:

Actively and passively bend and stretch the elbow.

#### Training:

Push-ups (optionally, against the wall)

#### Ankle:

#### Wrapping:

Fixation under the sole (foot bent).

Wrap over the ankle towards the leg.

Secure above the ankle. Optional: also wrap the heel.

#### Activation:

Actively and passively bend and stretch. "Chinese burn".

#### Training:

Kneeling lunge. Push knee forward.

### Vitality Flossing: „Compress to optimise“

This training method is based on compression-based techniques with a stretchy latex band, the Flossband. The Vitality Flossing concept focuses on dynamic muscle exercises, especially affecting joint mechanics and tissue smoothness.

#### SWATT method:

##### 1. Screening

Which movements or areas of the body are restricted?

##### 2. Wrapping

2.1 Pulling force: Apply the Flossband to the affected area, using 60 - 80 % tension (40 - 50 % on the opposite side).

Exceptions: Only 50 % tension for first aid & lymph wraps.

#### 2.2 Direction of pull:

Does the mobility of the affected part of the body improve if you move the surrounding tissue/skin outwards or inwards? This is the direction the Flossband should pull. Wrapping usually takes place from the extremity (distal) to the centre of the body (proximal).

#### 2.3 Fixation (first wrap):

Place the Flossband in the desired position and apply a complete wrap without tension with 100 % overlap.

#### 2.4 Wrapping:

Carry out the further wraps with 50 % overlap. You can also double-wrap the main pain area or use an X-shaped wrap.

#### 2.5 Securing wrap

The last wrap should be slightly looser so you can secure the end of the Flossband.

#### 2.6 Duration

Maximum 2 to 5 minutes per application. The greater the tension, the shorter the duration. You can repeat (even several times) after 2 to 4 minutes.

#### **3. Activation**

Active and passive movement of the wrapped parts of the body. Intensify the movement of the surrounding tissue by using a "Chinese burn" twisting motion.

#### **4. Training/Therapy**

Complementary closed-chain movements.

#### Wrist:

#### Wrapping:

Fixation on the palm of the hand.

Wrap twice over the back of the hand, then over the wrist.

Secure the band at the forearm.

#### Activation:

Actively and passively move the hand in all directions.

Support the wrist at a 90° angle.

#### Knee:

#### Wrapping:

Fixation below the knee.

Wrap over the knee towards the thigh.

Secure the ends of the band at the thigh.

#### Activation:

Actively and passively bend and stretch the knee. "Chinese burn".

#### Training:

Squats.

#### Shoulder:

#### Wrapping:

Fixation at the centre of the upper arm.

Wrap upwards until the shoulder joint is completely covered.

Secure the ends of the band.

#### Activation:

Move/rotate the shoulder joint whilst lying on your back.

#### Training:

Move shoulder back and forth whilst kneeling. Optional: push-ups.

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©2018 Ludwig Artzt GmbH. All rights reserved. Printed in Germany.