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Introduction

As a specialist for fitness and health products, Ludwig Artzt GmbH has been bridging the gap between sports science and the world of fitness for more than 35 years. Resistance devices constitute a particular focal point of our product range. Under the ARTZT vitality® umbrella brand, specialists and home users have access to a comprehensive selection of high-quality products and modern, efficient training concepts for professionals and recreational athletes.

Ongoing, close communication with experts is the pillar that supports the success of Ludwig Artzt GmbH. Our own training and continuing education institute is considered one of the key points of contact for qualified seminars on the topics of health and fitness. Disseminating technical knowledge backed up by sports science is the focus of our educational activities.

The elastic bands can be found exclusively in many exercise programs. Clinical research has proven that elastic resistive exercise increases/improves:

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Material and characteristics

- » High-quality, latex-free material (polyisoprene) from medical technology
- » Great grip thanks to special surface structure
- » Linear energy build-up that protects the joints and prevents injuries
- » Washable and easy to clean
- » Powder-free
- » Suitable for allergy sufferers
- » Available as a set with 2 exercise bands (medium and strong) and 1 Assist to easily attach the bands.

Enthaltene Widerstandsstufen/ contained resistance levels

mittel/medium

stark/strong

Factors that influence durability and safe exercising with the bands:

- » Avoid mechanical impacts (sharp-edged objects), as small holes or lateral tears can lead to the Exercise Band ripping.
- » Check the band prior to every exercise!
- » Use the ARTZT vitality Assist to attach the band. Do not use knots or plastic clamps.
- » Protect the Exercise Band from extreme temperatures and direct sunlight.
- » If possible, store the Exercise Band rolled up or hanging up.
- » Dirty Exercise Bands can be easily cleaned with soap and water and are ready for use again after they have dried.
- » Improper or careless use of the band may result in serious injury. Therefore only perform safe and painless exercises. Avoid exercises that could cause the band to snap back towards the head, or wear protective goggles during these exercises.
- » The Exercise Band is not a children's toy! Children should only exercise under supervision.
- » Avoid using bands that are too short. The ideal band length is 2.50 m. This prevents the buildup of high forces, uncontrolled movements and, above all, tearing of the bands.

- » The Exercise Band^{PLUS} can be used in swimming water. After use, the band should be rinsed with clear chlorine-free water and dried while flat.

Selection of the correct band


The appropriate colour / resistance of the band depends on your individual physical condition.

Choose a band that allows you to complete 15 repetitions.

Wrapping and connecting of the bands

Your elastic band should be securely attached to your hand or foot before use to avoid slippage and possible injuries. Double-wrapping the band may help to secure it to your hand or foot.

Attachement to the hands



Wrapping at one hand

Lay the band flat in your hand with the end towards your pinky finger. Wrap the long end of the band around the back of your hand. Repeat as needed. Grasp.

Wrapping with both hands

Begin with palms up and end of band between the thumb and finger. Rotate your palms inward, bringing the band around the back of your hands. Repeat as needed. Grasp.

Attachment to the feet

Foot Loop

Loop over the top of the foot and stabilize other ends with opposite foot.

Foot Wrap

Stand in the middle of the band. Wrap one end around the top of the foot.

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The main aims of a health-oriented preventive exercise program are increasing endurance, muscle mass and to improve bodyshape and posture. To match those goals, beginners should choose 6 to 8 exercises from the following program. Complete 1 to 3 sets of 15 to 20 repetitions with mild fatigue on the last set. At no time should you feel the strain of the exercise as very heavy; if so, change to an easier band colour. You may feel more comfortable using different levels of resistance bands for different body parts! Progress to the next colour of band or add sets or training units to your program when you feel comfortable and you are able to easily complete the 3 sets of 15 to 20 repetitions. Follow these guidelines when exercising with bands:

- » Perform motions slowly and under control. Maintain good overall body posture before, during and after each movement.
- » The band should be under slight tension before each movement to ensure the correct joint position.
- » The direction of the pull of the band influences both muscles and joints. Therefore it is important to be aware of the „force angle“ (angle between the band and the lever arm (arm, leg)). The maximum torque

(strength) is often noted at 90° (mid-range). Below 30° there is little significant effect.

- » Do not exercise while experiencing pain. If you feel pain stop exercising. If your pain should persist consult a health care provider.
- » All exercises are shown in the starting and final position. Closely examine the photos and imagine the movement. Always exercise BOTH sides of the body!

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In order to carry out pulling and rotation exercises effectively, the band can be attached to a stable anchor point (door handle, solid table/cupboard leg, etc.) using the ARTZT vitality Assist.

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1. Place the Assist around the anchor point.

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2. Pass the large opening of the Assist through the small opening and pull tight.

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3. Pass one end of the band through the large loop (photo 5) and pull it through to the middle of the band.

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Stand in a lunging position with a tense band behind. Grab the band with both hands with the elbows close to the head.

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Stretch elbows against the band's resistance and slowly return to the starting position.

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Stable parallel position. Close band position in front of the body. Elbows bent at about 90°.

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Raise arms slightly outwards against the band's resistance. Keep the lower arms parallel to the floor.

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Lunging position with a taut band behind. Grab the band's end at chest level. Keep the lower arms parallel to the floor.

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Push the hands forcefully forward until arms are extended.

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Grab the middle of the double-layered band with one hand behind your back at the level of your lumbar. The other hand holds the resulting band's loop next to the body.

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Push the band forwards forcefully. Switch hands.

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Stand with bent knees and slightly lean your upper body forwards. Stretch the band from the bottom-front.
Pull the ends of the band forcefully outwards and upwards.

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Stand with bent knees and slightly lean your upper body forwards. Stretch the band from the top front.
Pull the ends of the band forcefully from the bottom back.

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Stand with bent knees. The nearly stretched arms grip the band ends laterally in front of the body at shoulder height.

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Pull the band powerfully and diagonally across the other side with stretched arms.

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On knees, with a taut band from above. Grab the hanging band's ends from above.

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Pull the ends of the band/ elbows forcefully towards to the bottom-back.

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Straight standing position. The band is double layered, with the left foot in the resulting loop. Grab the band's end at hip level with the right hand.

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Spread the left leg counter to the band's resistance. Switch legs afterwards.

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Stand on the middle of the band. Grab the crossed band's end at hip level.

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Switching hands, pull outwards and upwards.

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VARIATION: Spread the opposite leg parallel to the pulling arm.

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Stand on the middle of the band. Grab the band's end at head level. Elbows bent.

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When standing up, stretch both arms above the head.

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©2019 Ludwig Artzt GmbH. All rights reserved. Printed in Germany.
Ludwig Artzt GmbH - Schiesheck 5 - 65599 Dornburg - Germany
©2019 Ludwig Artzt GmbH. Alle Rechte geschützt. Gedruckt in Deutschland.