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User manual Senso® Balance Pad **16 cm and 20 cm diameter**

Congratulations!

You are now the owner of a Senso® Balance Pad, a product “made in Germany”.
This sports tool was produced in compliance with valid EU norms. It is free of forbidden phthalates and heavy metals.

This tool is completing the Senso Balance product family. The soft and durable material and the knobs make a new level of sensorimotoric training possible. The pads are especially suitable for the introduction into walking and developing perception skills, but also for all sports and gyms. Because of the different colors available and the low depth they can be used as visual markers.

Warning

- Avoid storage in direct sunlight or near heat sources, because otherwise the SensoBalance Pad might deform
- Avoid contact with sharp and pointed objects (e.g. broken glass, pointed rocks in shoe soles) because they can damage the SensoBalance Pad. Do not throw the SensoBalance Pad!
- Use sport shoes, regular shoes or be barefoot while training on the Pads.
- Prior to any use, the SensoBalance Pad should be checked for damage because otherwise a safe handling is not guaranteed
- The SensoBalance Pad has to rest on a horizontal, stable and dry ground, to make sure, that the Pad cannot slide away during the exercises.

Equipment details:

The SensoBalance Pad is made of high-grade Rulon with a needle valve.

Size: 16 cm Ø, ca. 20 cm Ø

Care:

Use mild standard detergents and disinfectants. Dissolvers as well as acidic or corrosive cleansers can harm the material. The SensoBalance Pad is waterproof