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TOGU® Challenge Disc 2.0 Instruction Manual

Intended Usage: TOGU Challenge Disc can be used in athletic and fitness training to improve ones coordination and proprioception and in Physical Therapy for injury prevention and rehabilitation of ones ankles, knees, hips and spine.

Requirements: Constant supervision and instruction from a doctor, certified trainer or therapist is required when using the TOGU Challenge Disc for Physical Therapy.
Age: from 4 to 100 years, according to the level of fitness.
Weight: max. 120 kg body weight (approx. 270 lbs)
Skills: no skills necessary

Safety Instructions: One must execute all workouts on a compact, non slippery, dry and level surface must be clean and dry during all training. One can train with sneakers or barefoot. Do not wear socks, stockings or dress shoes. If necessary, one can use a chair, table, bar, or spotting person for assistance during fitness training. One can reduce the risk of injuries when falling by keeping the workout area clear of debris and if needed, place a soft surface around the training area. Children must only exercise under supervision of an adult.

Additional Instruction: The TOGU Challenge Disc must not be used for jumping, but only for smooth and easy movements.

Regulations/Maintenance: Before usage check the device for damages and make sure that all adjustable pieces are fixed.

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TOGU® CHALLENGE DISC 2.0

by

Software download & support - MFT Bodynetwork

mft-company.com

MFT Bodynetwork - Coordination/straining & Balance Boards

Contact Suppliers SOFTWARE DOWNLOAD & SUPPORT Deutsch

MFT PRODUCTS BODYTEAMWORK QUALITY ONLINE SHOP

MFT CHALLENGE DISC SOFTWARE

- WINDOWS AB 8.1, 10 Software Challenge Disc 2.0 Bluetooth + USB
- WINDOWS 7, VISTA, XP SP3 Software Challenge Disc USB
- APPLE MAC OSX, MACOS Software Challenge Disc 2.0 Bluetooth + USB

INSTRUCTIONS

- INSTRUCTION MANUAL (PDF) MFT Challenge Disc 2.0
- README MFT CHALLENGE DISC 2.0 (PDF) Notes and troubleshooting for the MFT Challenge Disc 2.0
- README MFT CHALLENGE DISC USB (PDF) Notes and troubleshooting MFT Challenge Disc USB

TROUBLESHOOTING AND MORE

- TROUBLESHOOTING MFT CHALLENGE DISC 2.0 (PDF) Diagnose and fix common issues – Bluetooth sensor does not connect with your tablet/PC
- HOW-TO VIDEO BATTERY CHANGE AND/OR SENSOR INSTALLATION Video instruction on how to change the battery and/or install the sensor of the MFT Challenge Disc

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Software-downloads:
www.challenge-disc.com
AppStore or GooglePlay

Bluetooth sensor activation

Please briefly move / tap disc
Start software

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Calibration

When starting the software programme, the Challenge Disc must be lying on the floor. You must not stand on the disc at this stage! The disc's position should not be changed during training, otherwise it may be necessary to recalibrate the disc, which would include restarting the programme as under point 1. Every restart automatically leads to a recalibration. The disc is positioned correctly if the USB cable is pointing towards the laptop.

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