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
## Introduction

Posture and stabilisation are terms we encounter with increasing frequency in relation to muscle and joint complaints. But what does “posture” actually mean? According to the sports medicine definition, it refers to maintaining an upright, functional and economical body position that counteracts the effects of gravity. Since we are always working against gravity, ongoing muscle activity is required in order to stabilise the joints and extremities in continuously changing positions without unnecessary fluctuations / avoidance movements or excessive muscle activity (tension).

These ongoing corrections, which are required in order to stabilise the upright position, are calculated and controlled by the central nervous system with the help of synergistic muscle activation, i.e. the simultaneous activation of muscles that have opposing effects on a joint such as the arm flexors and arm extensors. When this control process is working smoothly, the body is in “balance” – it is stable and can execute purposeful movements.

Monotonous positions such as sitting and standing or continuously recurring movement patterns frequently dominate our everyday lives. The lack of breaks or compensating movements cause the “stabilisers”, which are

mainly controlled by the brain, to become fatigued and the global muscles (“movers”) to assume their function. But they are not intended for this purpose over the long term and respond with tension. This results in muscular imbalance with altered tensile and compressive loads on the joints, which can lead to long-term postural deformities and injuries. In addition, the ability to activate the deep stabilisers is lost.

The ARTZT vitality® Bioswing® is a swing stick with excellent swing characteristics (quality of the oscillations) which facilitates verification of stabilisation efficiency as well as individual training with a postural effect (stabilisation of posture and movement). With regulated impulses of great stimulation density, the exercises with the Bioswing® systematically promote the interplay of stabilising musculature and thereby establish the foundation for purposeful movements and a stable body posture.

Tense muscles loosen up again and remain relaxed even under strain. Your movements become more precise and can be carried out without an excessive expenditure of energy. In short – your pain can go away and you experience newfound success in your training!


## Areas of application

The Bioswing® is used in neuro-orthopaedic rehabilitation, for example in case of chronic muscular imbalances or after injuries in the shoulder and spine region.

Athletes who are subjected to one-sided strain (e.g.: tennis, spear throwing, golf) use the ARTZT vitality® Bioswing® to strengthen their muscles and stabilise their posture.

On a day-to-day basis, people involved in occupations that mainly involve sitting down use the Bioswing® to conduct effective full-body training through the neuro-muscular activation of the body's key regions (the spine, shoulder girdle and pectoral girdle).

The neuro-physiological effect is reflected by the fact that incoming impulses of regulated stimulation density cause the simultaneous activation of muscles that have opposing effects on a joint (synergistic muscle activation). This effect is of crucial importance for any kind of stabilisation demands on the locomotor system. Therefore motor function and vertical stability can be systematically optimised with the Bioswing®.

## Materials and characteristics

150 cm long and a mere 627 g in weight, the ARTZT vitality® Bioswing® is made of high-tech spring steel which was especially tested and is also used in the automobile industry. This material does not break and is oscillation elastic. The special feature of the Bioswing®: dampening elements on the handle mean the swing stick responds to movements smoothly. This reduces the direct load on the wrist and therefore on the basal thumb joint. Uncontrolled movement of the swing stick and impacts on the musculature are prevented. The Bioswing® is also equipped with individually adjustable frequency regulators so that the user can continuously adjust the oscillations.

The Bioswing® from ARTZT vitality® supports the continuous adjustment of the exposure intensity from 3.2 to 4.4 hertz (Hz). In practice, this means: the further out the frequency regulators are positioned, the easier the exercises are to complete. This means one and the same swing stick offers effective training for anyone from beginners to ambitious athletes. The Bioswing® was developed in the course of long-term cooperation between Haider Bioswing and Univ. MUDr. Med. Eug. Rasev (professor and chair for rehabilitation and physiotherapy at Charles University in Prague and founder of "postural therapy") based on the scientific insights of Prof. Vladimír Janda, MD, DSc.

## Handling

Different parts of the body can be stimulated by choosing the swing stick position with a one-handed or two-handed grip and switching between defined body postures and oscillation planes.


## Oscillation planes

The Bioswing® can be swung in two planes.

- » back and forth (ALWAYS as an extension of the forearm)

- » up and down OR left to right (ALWAYS at a 90-degree angle to the forearms)

## Grip types

The Bioswing® may be held with one or both hands. In general, you should never grip the swing stick with maximum force (clamping down on the hand grip) but as loosely as possible while holding the Bioswing® securely (it should not slide through your hand).

When gripping the swing stick with one or both hands so it is parallel to the floor and when gripping it with one hand and the swing stick in the vertical position, the wrists should always be angled slightly towards the back of the arms. With the two-handed grip, both hands meet in the middle (see illustration). The one-handed grip is applied in the middle of the hand grip.

When gripping the swing stick with both hands so it is perpendicular to the floor, the hands can either be loosely positioned one over the other (see illustration) or folded (prayer position).

In addition to the factors already mentioned above, the durability of and safe training with the ARTZT vitality® Bioswing® mainly depend on proper use:

- » The Bioswing® should be stored hanging up. Alternatively you can place the swing stick flat on the floor. When doing so, the frequency regulators must be moved to a centred position so that the swing stick does not sag. Over the long term, sagging can cause the swing stick to bend slightly which has a negative impact on optimum oscillation quality. The swing stick should only be leaned upright against a wall for short periods of time.
- » Improper or careless use of the Bioswing® may lead to injuries. As a result, only safe and pain-free exercises should be completed.
- » The Bioswing® is not intended as a toy for children! Children should be supervised when completing exercises. In general, effective training is only recommended with children aged approximately 12 years and up.
- » Please observe the amount of space required to train with the Bioswing®: approximately 4 m<sup>2</sup> per participant.
- » When completing exercises with the swing stick in the vertical position, you should keep the swing stick tilted slightly and / or reduce the oscillation range so that the ends of the swing stick do not strike the head or legs.

### Body posture

In order to assure the effectiveness of the following exercises, understanding and assuming the best possible body posture is essential. This relieves stress on the passive structures (bones, tendons, ligaments) while assuring optimum training of the targeted musculature. This ideal body posture (while standing up) is also defined by conscious alignment of the head, shoulder girdle, spine, pelvis and feet together against gravity. Please pay particular attention to assuming this posture!

### In general:

- » Stand on both legs with your feet shoulder-width apart
- » The back of the head, shoulder girdle and buttocks form a line..

### Head position:

- » Looking ahead and down.
- » Push the back of the head backwards without lifting the chin (double chin).

### Upper body:

- » Slightly raise the breastbone forward (“show off your medal”).
- » Allow the shoulder blades to drop back and down (towards your buttocks).
- » Fluent breathing (mainly abdominal and flank breathing).

### Legs:

- » Hip and knee joint slightly bent. Upright torso, leaning ahead SLIGHTLY.
- » Slight exterior rotation of the legs.
- » Straight leg axis (no O or X leg position).

### Feet:

- » Slight exterior rotation of the feet (10-15 degrees).
- » Balanced load on the foot (forefoot and hindfoot).
- » Barefoot if possible.

### **Training scope**

Precisely adapting the duration of training and allowing sufficient breaks (loosening / stretching) are important for effective stabilisation training.

For this purpose we differentiate between a training unit, exercise repeat and training cycle. A training unit consists of the uninterrupted training time. This should be between 10 and 15 seconds. Counting to 10 in your head is a good indicator.

An exercise repeat consists of repeating the training unit once. The optimum exercise repeat count is between 5 and 10.

A training cycle consists of completing all exercises once with the selected repeat count (5-10). 1 or 2 training cycles are recommended. This corresponds to a maximum training period of 20 to 30 minutes.

### **General structure of a training unit**

Ideally you should warm up before training with the Bio-swing® and prepare your body for the stabilisation training that will follow (max. 15 minutes). A general warm-up may consist of cycling, jogging (running in place is also possible) or jumping rope.

If you feel very tense or strained, you should definitely do some stretches (neck, chest, hips).

The specific warm-up or purposeful preparation may consist of the following exercises:

1. Describe circles with your open hands (fingers spread).
2. Stand on one leg: describe circles with your ankle (toes pulled back) and lower leg. Ensure your hip position remains stable!
3. Describe backwards circles with your shoulders and arms. The thumb leads the movement while the elbow is bent over the course of the entire circle. The palm is open and the torso follows the arm movement slightly (rotation).
4. Back of the head stretch: loosely clasp your hands behind your head at the level of the hairline. The elbows face forward without pulling on the head. Now push against your clasped hands with your head (tuck in the chin / double chin position). Hold briefly and then return to the initial position. The movement is NOT a tilt but a "horizontal push".

5. Push your head into the palms as described under step 4. and go into a squat with your upper body straight and leaning ahead slightly (keep your buttocks back as though you were about to sit on a chair!). Only go so low that the entire soles of your feet stay on the floor and your back remains straight.

The main component with the Bioswing® exercises follows the warm-up. It is important to start slowly, with only a few training units (3-4). In doing so, observe the following information:

- » You should take an active break (5-10 seconds) by shaking out and relaxing your arms and legs after each training unit.
- » Maintain a stable head, shoulder and pectoral girdle position (i.e. nearly motionless) and be sure to keep your shoulder position low!
- » If you initially experience problems completing the entire exercise, break the exercise down into individual training steps (e.g. below shoulder height, slightly above shoulder height, over your head). Gradually combine the steps you are able to complete smoothly and without avoidance movements until a flowing transition develops.
- » If you initially find it difficult to maintain a stable pelvis position, complete the exercises in a slight walking position ("anchoring" the pelvis). As you advance, you can make the exercises more difficult by changing the position to the parallel stance (narrow / wide) and the one-legged stance or by means of additional visual stimuli (closing your eyes).
- » First build up the number of repeats.
- » Only then should you move the frequency regulators further in step by step whenever you can complete 5 to 6 repeats of an exercise readily and with confidence.
- » Training should be pain-free. Take a break if an exercise becomes painful. You should consult a physician if pain occurs repeatedly.
- » Study each exercise in detail and visualise the process. ALWAYS exercise both the right and left side of the body even when the illustration only shows one side.

### **Discontinuation criteria and signs of excessive strain**

If you experience any of the following "symptoms", you should discontinue training or only begin training on the advice of a doctor or therapist:

- » Ongoing pain.
- » Increased concurrent movements of the hand, shoulder girdle or pectoral girdle.
- » Muscle cramps.
- » Inability to keep the Bioswing® oscillating.
- » Reaching the maximum training duration.
- » Painful movements of the cervical spine.
- » Limbic dysfunction (emotional stress).

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Hold the Bioswing® in front of the body with one hand above the other (perpendicular to the floor). The elbows are bent.

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Now make the swing stick oscillate by moving your hands forward and backward forcefully and keep it oscillating while turning the upper body slightly to the left and then past the centre and to the right (approximately 5-10 degrees). The head follows the movement of the upper body.

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**IMPORTANT:** Move the upper body only while keeping the pelvis stable. Be sure to maintain a low oscillation amplitude (20-30 cm) during this exercise to avoid collisions between the swing stick and your body. Also switch your hands after each training unit (upper hand down – lower hand up).

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
**IMPORTANT:** Do not pull up your shoulders at any point during the movement and always keep your arms slightly bent.

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Hold the Bioswing® at the level of your navel with both hands (parallel to the floor). The thumbs face each other and touch lightly. The elbows are bent.

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Now make the swing stick oscillate by moving your hands forward and backward forcefully and keep it oscillating. Describe a semicircular movement with your arms until they are above your head (elbows approximately at ear height) and then lower the arms back down to the level of your navel.

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Go into a squad (approximately 110 degrees) and bend the upper body forward slightly. Hold the Bioswing® just above the knees with both hands (parallel to the floor). The thumbs face each other and touch lightly. The elbows are bent and slightly more of your weight is shifted to the forefoot.

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Now make the swing stick oscillate by moving your hands forward and backward forcefully and keep it oscillating. The pelvis remains stable.

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Hold the Bioswing® with one hand (from above) at navel height (diagonally in front of the centre of your body). The back of your hand points up and out, the elbow is slightly bent. The free hand faces forward with the palm and the fingers spread.

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Now make the swing stick oscillate by moving your hand forward and backward forcefully and keep it oscillating while moving the swing stick out, up and back in a semicircular movement.

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
In the final position, the back of the hand is facing down while the hand and shoulder are in one line (with the elbow slightly bent).

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**IMPORTANT:** Keep your shoulders low and the elbows bent. The upper body remains stable and does not sway. If you find this exercise easy, you can try moving your arms forward and backward.

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Hold the Bioswing® at chest height with both hands (parallel to the floor). The thumbs face each other and touch lightly. The elbows are bent.

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Now make the swing stick oscillate by moving your hands up and down forcefully and keep it oscillating.

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
**IMPORTANT:** The arm and leg positions remain unchanged during the exercise. After each training unit, switch the supporting and non-supporting leg.

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One-legged stance with the hip joint bent approximately 30 degrees and slight exterior rotation of the thigh. The lower leg “hangs” down vertically and the toes are pulled back. Hold the Bioswing® at the level of your navel with both hands (parallel to the floor). The thumbs face each other and touch lightly. The elbows are bent.

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Now make the swing stick oscillate by moving your hands forward and backward forcefully and keep it oscillating.

# UPGRADE YOUR VITALITY!

**ARTZT**  
vitality

## Produktprogramm ARTZT vitality® | Product range ARTZT vitality® Medical Fitness - Made in Germany

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
- **Flowin® Personal Training**
- **Gymstick™ Health**
- **BOSU® Balance Training**

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ARTZT vitality® is a registered trademark of Ludwig Artzt GmbH. Unauthorized use is strictly prohibited.  
© 2009 Ludwig Artzt GmbH. All rights reserved. Printed in Germany.

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© 2009 Ludwig Artzt GmbH. Alle Rechte geschützt. Gedruckt in Deutschland.